


Interpretation guide

ImmunoCAP® Rapid offers easy-to-interpret results

Allergen-specific results reported as Class 1, 2, or 3

ImmunoCAP Rapid Inhalant Profile 1 reports



Results are reported in classes (see table below).


Summarized result interpretation

ImmunoCAP Rapid Inhalant Profile 1 result	Level of allergen-specific IgE antibody (approximate kU _A /L)	Degree of sensitization to specific allergen
Class 1	<1	Absent or low
Class 2	≥1 and <15	Moderate to high
Class 3	≥15	Very high

With ImmunoCAP Rapid results in hand, the clinician can select appropriate medications and target exposure reduction measures based on the individual's allergic sensitization.

ImmunoCAP[®] Rapid results guide targeted disease management

Management options according to the allergic response for inhalant allergens*



*Results must always be interpreted based on clinical signs and symptoms.