



## Quality Club Specific IgE General Monthly Report


Test Period: August 2012  
Sample Lot No.: AWXB7

### Results from all labs

N (except outliers)	461
Median	7,98 kU <sub>A</sub> /l
Mean	8 kU <sub>A</sub> /l
SD	0,89 kU <sub>A</sub> /l
CV	11,2 %
Mean - 2 SD	6,22 kU <sub>A</sub> /l
Mean + 2 SD	9,79 kU <sub>A</sub> /l

### D. farinae (d2)

No of participating labs: 463




### Results from all labs

N (except outliers)	455
Median	1,12 kU <sub>A</sub> /l
Mean	1,13 kU <sub>A</sub> /l
SD	0,14 kU <sub>A</sub> /l
CV	12,3 %
Mean - 2 SD	0,85 kU <sub>A</sub> /l
Mean + 2 SD	1,41 kU <sub>A</sub> /l

### Egg white (f1)

No of participating labs: 464



### Results from all labs

N (except outliers)	415
Median	30,9 kU <sub>A</sub> /l
Mean	31 kU <sub>A</sub> /l
SD	3,21 kU <sub>A</sub> /l
CV	10,4 %
Mean - 2 SD	24,6 kU <sub>A</sub> /l
Mean + 2 SD	37,4 kU <sub>A</sub> /l

### Betula verrucosa, Silver birch (t3)

No of participating labs: 419

